

How to Download a Copy of the Antech Order Codes with Ext/Profile ID's, to add to your Treatment List in AVImark (to use the Antech Lab Requisition Form)

AVImark will automatically enter the Antech Ext/Profile ID when the treatment is defined, but it is the Clinics responsibility to enter the Antech Order Code on the Advanced tab of the Change Treatment window within the Treatment List of AVImark. See screen shot below for an example of where to enter the above information.

The screenshot shows the 'Change Treatment' window with the 'Advanced' tab selected. The 'Antech Ext/Profile ID' field is set to '921' and the 'Antech Order Code' field is set to 'LD50'. Red arrows point to these two fields. Other fields include 'Vaccine' (none), 'Remind As' (none), 'VetConnect ID' (NONE), and 'Weight Surcharge' (0%). The 'Instrument Types' list includes 'None' (checked), 'Abaxis HMT', 'Abaxis Vet Scan', 'CDC Hemavet', 'Hemagen', and 'Heska ABC'. The 'Follow-up' section has 'Days' set to 0 and 'Follow-up Critical' unchecked. The 'Company' field is empty. At the bottom, there are 'Next', 'Prior', 'OK', and 'Cancel' buttons.

Both AVImark and Zoasis have a list of the most commonly used Antech Order Codes, with their Antech Ext/Profile ID available on their websites. There is a separate list for Antech East clients and Antech West clients. This list will be updated on a bi-annual basis.

1) You need to recognize if your clinic is an Antech East or Antech West customer. Many of the codes are different depending on which lab your hospital uses. If you are unsure which list applies to your hospital, you can quickly find out by looking on your Antech Test Request Forms.

Antech West clients will have a 1W near top right of form:

The form header for Antech West clients includes the Antech Diagnostics logo, 'LAB USE ONLY', and 'PLEASE LABEL ALL TUBES WITH PT. NAME & ACCT. NO.'. On the right side, there is a 'CRITICAL REORDER FORM' checkbox and a '1W' in a circle, with a red arrow pointing to it.

Antech East clients will have a 1E near top left side of form:

The form header for Antech East clients includes the Antech Diagnostics logo, 'LAB USE ONLY', and 'PLEASE LABEL ALL TUBES WITH PT. NAME & ACCT. NO.'. On the left side, there is a '1E' in a circle, with a red arrow pointing to it. On the right side, there is a 'SPECIMEN REQUIREMENTS' table.

SPECIMEN REQUIREMENTS			
S Serum	C Culturette	R Red (non SST)	SU Unopened Serum
L Lavender	SL Side	B Blue	BCB Blood Culture Bottle
U Urine	F Feces	G Green	FSS Frozen Serum Separator

2) Download the appropriate list of Antech Order Codes (East or West) from either AVImark.net OR Zoasis.com (if using Zoasis to download list, you will need to be logged into Zoasis as your hospital user in order to access it.)
www.avimark.net - Go to Technical Support > Downloadable Documents > Antech Order Codes.
www.zoasis.com - Login, go to Laboratory > Test Guide > Antech Order Codes w Ext/Profile ID's.

Both sites have the Antech Order Code lists in PDF format. You will need Adobe Acrobat Reader to access these lists; if you do not have it you can download a free version from Adobe at <http://www.adobe.com/products/acrobat/readstep2.html>. You can either save the PDF file to your computer or you can print it (by clicking the Save or Print options from within Adobe.)

3) Once you have the list of Antech Order Codes and their corresponding Antech Ext/Profile ID's, you can quickly figure out what Order Codes your hospital most commonly uses by printing a copy of your Top Tests Business Report from Zoasis.com.
 -Log into Zoasis.com and go to Laboratory > Reports. Locate the Business Report called Top Test Report. You can print the report by clicking on Actions > Print.

The screenshot shows the Zoasis.com Laboratory Reports interface. The 'LABORATORY' menu is highlighted with a red arrow. The 'REPORTS' menu is also highlighted with a red arrow. The 'Top Tests' report is displayed for the month of July 2007. The table below shows the top tests with their respective metrics.

Order Code	Description	Report Month	Previous Month	Change (%)	YTD		
405	Ova & Parasite/Giardia (ELISA)	382	419	-8	2374	0	0
T615	Occult Heartworm Antigen	181	231	-21	1211	129	838
5013	Geriatric Senior Profile	128	117	9	926	750	23
T835	Urine Microalbumin (Canine)	124	126	-1	858	619	38
85891	Urine Microalbumin (Feline)	120	96	25	767	657	16
80030	Geriatric Profile (No UA)	85	102	-16	619	618	0
3004	Profile Urinalysis	36	35	2	230	164	40
T495	T4	35	24	45	194	219	-11
CYTO	Cytology	34	41	-17	330	187	76
9978	Basic Wellness Screen	33	32	3	196	128	53
T760	Urinalysis	32	32	0	225	156	44
M130	Culture (Urine)	31	35	-11	148	0	0
99780	Basic Wellness Screen (No Urine)	27	43	-37	217	191	13
85424	Canine Wellness Profile	26	22	18	138	149	-7
FBX	Biopsy	17	13	30	114	143	-20
9523	General Senior Profile	14	8	75	51	94	-45
SA440	Hyperthyroid Panel	14	18	-22	106	84	26
T330	Complete Blood Count	14	12	16	65	66	-1
9816	Free T4 (Add On)	14	10	40	83	27	207
SA025	Vet-Screen	13	11	18	67	23	191

If you are still unsure about which Antech Code you want to order, contact your local Antech Customer Service team for help figuring out which code is for the test you are looking for.
 Antech East Customer Service: 800-872-1001
 Antech West Customer Service: 800-745-4725
 Antech Test Express: 800-EXP-TEST